

Frequency Electrical Stimulation by Sedatelec

AGISTIM DUO



AGISTIM DUO

Professional, multi-functional by design



Nogier Frequencies
Frequencies from
1-100 Hz



Timer
treatment modes



4 channels,
up to 16 needles

To enhance the effects of acupuncture treatments, effectively relieve pain or reduce localised unsightly lipodystrophies, **AGISTIM DUO** is your immediate choice.

Efficacy through SCFO* innovation

An ingenious program of permanent impulse frequency fluctuation around the chosen frequency greatly reduces the risk of physiological habituation to stimulation.

**Statistically Controlled Frequency Oscillation: a Sedatelec exclusively patented process*

Simplicity

- 2 independent simulators in the same instrument (Duo function)
- Individual keys to set time, frequency and intensity
- Multi-lingual control screen (EN, FR, ES, PT, DE, IT, NL)
- 4 independent channels
- 2 to 16 colour connectors
- Stimulation through **needles** and adhesive skin **electrodes**
- Timer, up to 60 minutes

Maximum safety

- Total insulation from the main supply
- Safety function prevents treatment beginning without resetting the intensity to zero.
- Stop alert if a frequency parameter or emission mode is changed during treatment
- 4 electrically insulated channels



Thanks to their original design, the needle connectors ensure perfect contact with the needle, with at the same time guarantee of complete insulation of the skin.

During treatment, the patient can adjust the intensity of each channel by himself as no incorrect manoeuvres are possible.

The Sedatelec signature: Nogier frequencies

The 7 pre-programmed Nogier frequencies (from A to G) can be selected individually or in scanning mode for auriculotherapy or acupuncture stimulation.



AGISTIM DUO

Excellence in its field

Multiple functions, design, simplicity, safety and efficiency are the features of the AGISTIM DUO which you and your patients will appreciate on a daily basis.

**This confirms
your professional
approach and
the quality of your
treatment.**



Ask for the «For more information» documentation about AGISTIM DUO from your supplier or by writing to sedatelec@sedatelec.com

LAGISTIM DUO is not indicated for use in myostimulation rehabilitation. No FDA approved (USA).

Your distributor



Sedatelec • Chemin des Mûriers • F-69540 Igny • Lyon • France • Tél. +33 (0) 472 663 322
Fax + 33 (0) 478 508 903 • sedatelec@sedatelec.com • www.sedatelec.com

1 - 9,9 Hz

10 -

1 - 9,9 Hz

10 -