

premio 32 evo

Laser A-G program

NEW!

Effective,
Fast &
Deep
Stimulation



premio 32 evo

Laser A-G program



The Sedatelec Signature: The Nogier Frequencies

The 7 frequencies can be selected individually, or associated to develop the potential desired effects with the 4 available programmes:

— ANATALGY: E and G

Reduces the transmission of pain through the central nervous system. Its local action on the inflammatory process makes the patient feel better.

— TISSUE REGENERATION: A, B and F

Accelerates tissue reconstruction, restores cell functions, accelerates the scarring and regeneration processes both on the surface and within.

— MUSCLE RELAXATION: C, D and G


Eases release of tension and contractions.

— NOGIER FREQUENCY SCANNING

Stimulates the pathological zones regardless of the origin (pain or functional disorders, etc.) using successive emissions of all the Nogier frequencies.

Characteristics


Wavelength




Infrared **emission at a wavelength of 905 nm** passes through the cutaneous barrier to reach deep cells for safe penetration, as it is less absorbed by the skin's water, haemoglobin and melanin.

Peak power

This pulse mode, with a peak power of **100 Watts**, gives the laser beam **the power needed to penetrate to deep cells** and therefore to impart information that can be used for the effective local mobilisation of the organism's resources.



Superpulsed emission



Short pulse bursts (<100ns) have the tempo of the **Nogier frequencies**.

They impart **5 Joules** of energy per minute for 100 Watts of power with no destructive thermal effect.

- ✓ **CORDLESS:** handheld device without restrictions
- ✓ **EASY TO USE:** direct access to 4 programmes
- ✓ **PRATICAL:** the frequencies used are displayed, reference beep (30 sec)
- ✓ **PRECISE:** lights the zone being treated
- ✓ **ADAPTED TO YOUR PRATICE:** for local stimulation (ear, face and body) 



40 years of cold laser experience

Help your patients
with the latest in
Bio-Stimulation Laser technology

Case including:

- The PREMIO 32 EVO
- 1 power reduction nozzle
- 1 emission tester
- 1 pair of protective goggles
- 1 charger
- 1 user guide



*Prior to any use, carefully read
the user guide*



*Laser device, IEC 60825-1 class 3B
Invisible laser radiation.*

Exposure to the beam can harm the eyes.

Caution: Federal law restricts this device to sale by or on the order of a physician.



Sedatelec guarantees its products for 2 years and provides after-sales services from its premises in Lyon for up to 10 years after sales have been discontinued.

distributor's stamp

CE 0123
Made in France

FCP-SED-PREMIO32EVO-EN-2026-06 / © Sedatelec



Chemin des Mûriers • 69540 Irigny (Lyon) • France
Tel: +33 (0) 472 663 322 / Fax : +33 (0) 478 508 903
contact@sedatelec.com • www.sedatelec.com